## SMILE FALL

### User manual

NE41 19004-00 v2.0

Information in this user manual is subject to change without notice.

NEAT Electronics AB reserves the right to change or improve their products and to make changes to the content without obligation to notify any person or organization of such changes or improvements.

NEAT Electronics AB is not responsible for any loss of data, income or any consequential damage whatsoever caused.

For more information, details and descriptions, visit our web site:

[www.neat-group.com](http://www.neat-group.com)

© 2019 Legrand  
All rights reserved.

Document number: NE41 19004-00 v2.0  
Revision date: 2020-10-21

NEAT Electronics AB  
Varuvägen 2  
SE-246 42 Löddeköpinge  
Sweden

Phone: +46 (0)46 707065  
Fax: +46 (0)46 707087  
[www.neat-group.com](http://www.neat-group.com)  
[infosweden@neat-group.com](mailto:infosweden@neat-group.com)

#### Declaration of Conformity

Hereby NEAT Electronics AB declares that the radio equipment SMILE FALL is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address:

[www.neat-group.com/es/en/downloads](http://www.neat-group.com/es/en/downloads)

## 1 About SMILE FALL

The SMILE FALL is intended for users with low activity level and a risk of falling. At a high user activity level the battery life will be shorter than specified. The SMILE FALL is to be worn by the user on a wristband.

The SMILE FALL has the same general functionality as a SMILE, see NE41 14002-02.

### 1.1 Fall detection

SMILE FALL is developed to detect hard falls in situations where the user cannot press the button after the fall, including:

- slipping and falling backwards.
- falling due to loss of consciousness.
- tripping on an object and falling forward.

In case of a fall, SMILE FALL will detect it and send a fall alarm.

### 1.2 Other alarms

Falls where the impact is softened because the person manages to react and partly catch himself with his hands or knees, or because an object (e.g. wall or chair) is partly hindering the fall, will not be detected by the SMILE FALL. However, the

user can always raise an alarm manually with a button press.

## 2 Getting started

To activate the device, press the button and hold until the device starts blinking green. No fall alarms will be transmitted before activation has been done. If the device has already been activated a regular alarm will be sent when the button is pressed, indicated by red LEDs.

To pair the SMILE FALL with a radio receiver follow the instructions of the receiver.

A SMILE Programmer is required to customize the configuration of the SMILE FALL. Please contact your distributor.\*

## 3 Supervision

The receiver can be configured to supervise if SMILE FALL is in or out of range of the receiver. The SMILE FALL sends an auto test alarm for this purpose. SMILE FALL will also send low battery alarm as well as indicate battery status when sending any alarm.

\* Programmer not available in all markets

## 4 LED Indications

There are four LEDs located around the alarm button.

| LED indication                               | Event  |
|--|--|
| 1 red blink (during transmission, until ACK) | SMILE FALL transmits an alarm (triggered either from a fall or a push on the alarm button) |
| 1 green blink (2 sec)                        | The alarm is acknowledged (ACK)  |
| Flashing red or green (during transmission)  | Low battery, replace the unit  |

## 5 Technical data

|  |                                |
|--|--------------------------------|
| Measure                                | 34 x 42 x 15 mm                |
| Weight (incl. battery)                 | 14 g                           |
| Frequency <sub>RF</sub> - EU           | 869.2 MHz (Social alarms)      |
| Frequency <sub>RF</sub> - Non EU*      | 866, 868, 906, 916, 921 MHz    |
| Transmit power <sub>max</sub> - EU     | 869.2 MHz 10 mW (EIRP)         |
| Transmit power <sub>max</sub> - Non EU | According to local regulations |
| Battery life**                         | 1 year (on average)            |
| IP Class                               | IP67                           |
| Temperature range                      | +5 - +55 °C                    |

\* According to local regulation for social alarms.

\*\* At 4 alarms/day.

## 6 Important

- Always test the system per instructions prior to use.
- The product may not be suitable for all persons.
- Check device regularly and replace when necessary.
- The SMILE FALL is not intended as a life supporting device.
- Keep away from dust, moist and dirt.
- Do not warm up the device or use it near fire.
- For repairs, contact a NEAT dealer.
- Clean the device with a soft cloth, dampened slightly with mild soapy water.
- Do not clean the device with harsh chemicals, solvents or other corrosive substances.
- Do not open, no user serviceable parts inside.
- This product contains a lithium battery, swallowing battery may cause serious injury or death.
- Keep battery away from children.

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[www.neat-group.com/de](http://www.neat-group.com/de)

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NEAT GmbH  
Carl-Zeiss-Ring 12  
85737 Ismaning

Tel: +49 (0)89 519966 6-0  
Fax: +49 (0)89 519966 6-25

[www.neat-group.com/de](http://www.neat-group.com/de)  
[infogermany@neat-group.com](mailto:infogermany@neat-group.com)

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[www.neat-group.com/de/downloads/dokumente/](http://www.neat-group.com/de/downloads/dokumente/)

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